



Electric Examination
Table
Elevate Comfort
with Our Premium
Hospital Bed



SCAN FOR VIDEO
www.procaremedfurn.com

	PRO076 (Electrical examination couch)	
No. of Motor	1 motor	
Total Length	72 Inch	
Total Width*1	24 Inch	
Total Height	17 Inch to 31 Inch	
Mattress Deck	35 cm to 70 cm	
Safe Working Load*3	130 kg	

CE

*1 When the side-rail is raised.
*2 Indicates the height from floor to the mattress deck surface.
*3 Maximum load with which the bed can be operated. (The total load of user's weight and incidental items including the mattress and options)

Procare

Address: Gate no 197, dehu alandi road,
jadhavwadi, chikali Pune,
Maharashtra-411062
Mob no: +91 8956584884,8956114459
Email id: Info@procarestore.in
GST NO: 27ABCFP8176F1ZY



Accompanying the Recovery Journey

Our Electric Examination Couch is designed to offer maximum comfort and convenience for both patients and healthcare providers. Ideal for clinics, hospitals, and examination rooms, this couch is equipped with advanced features to facilitate smooth and efficient medical examinations.



Height Adjustment



Battery Backup




IV Rod Provision

Create a Safer Environment

KEY FEATURE

Electric Height Adjustment: Easily adjustable height to accommodate various examination needs, ensuring comfort for both the patient and the medical professional. It also comes with a optional battery backup feature.



Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare