



Create a Safer Environment

Accompanying the Recovery Journey

EXAMINATION COUCH 1 DRAWER 1 CABINET

SKU: PRO075

www.procaremedfurn.com

- Length: 1830 mm
- Width: 610 mm
- Height: 810 mm

KEY FEATURE

Premium Construction:

The couch is made from premium CRCA steel that offers the best durability, and flexibility to the Couch to use in daily medical examination procedures.

Pretreated and Epoxy Powder Coated:

The examination couch is protected with an epoxy powder coating that resists corrosion and provides a shiny and easy-to-clean surface which helps to maintain hygiene.

Multiple Storage Option:

This examination couch comes with a sufficient storage option which allows for storing emergency medicines and required equipment. This Couch comes with One drawer and One Cabinet.

Two-inch mattress:

This examination Couch comes with a 2-inch premium mattress to make this couch more comfortable for the patient

Foot Step Stool:

A Single foot step stool accessory is provided with the couch, Which helps to access the couch easily.

Premium Quality Cushion



Step Stool →

Procare

Address: Gate no 197, Dehu Alandi road,
Jadhavwadi, Chikali Pune,
Maharashtra-411062
Mob no: +91 8377832832
Web Site : www.procaremedfurn.com
GST NO: 27ABCFP8176F1ZY

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare

CE