



FIVE-FUNCTION  
MANUAL ICU BED  
Elevate Comfort  
with Our Premium  
Hospital Bed



## Accompanying the Recovery Journey


Recovering patients require a nurturing environment where they can recuperate with peace of mind. Studies have shown that a substantial number of adverse medical events involve falls. Placing particular emphasis on safety, the Qualitas Nx bed is designed to reduce the risk of incidents and provide superior care.



Ergonomic Side-rail Integrated Panel



Easy Grip Handle for Transportation



Full Length Side-rail

Create a Safer Environment

SCAN FOR VIDEO  
[www.procaremedfurn.com](http://www.procaremedfurn.com)

	PRO004 (Manual Hospital Bed)	
Total Length	228.6 cm	
Total Width*1	106 cm	
Total Height	73.66 cm to 104.14 cm	
Mattress Deck	35 cm to 70 cm	
Height*2	90.17 cm x 205.74 cm	
Safe Working Load*3	230 kg	
Caster	12.5 cm dia. Single-wheel	
Head of Bed Angle	0° to 75°	
Knee Angle	0° to 45°	
Trendelenburg	12°	
Reverse	12°	—
Trendelenburg		



\*1 When the side-rail is raised.  
\*2 Indicates the height from floor to the mattress deck surface.  
\*3 Maximum load with which the bed can be operated. (The total load of user's weight and incidental items including the mattress and options)

# Procare

Address: Gate no 197, dehu alandi road,  
jadhavwadi, chikali Pune,  
Maharashtra-411062

Mob no: +91 8956584884,8956114459

Email id: [Info@procarestore.in](mailto:Info@procarestore.in)

GST NO: 27ABCFP8176F1ZY


## Fall Prevention

Accidents do occur, even in hospitals. Studies have shown that a substantial number of hospital-acquired conditions involve falls. Hence, PRO is provided with integrated features reducing the patient fall risk.

Ratio of Fall Incidents in all  
Medical Adverse Events


%19  
Falling down  
2%  
Falling over

Other  
Incidents



Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Infection Control



Removable Mattress Deck and Smooth Surface for Easy Cleaning

The mattress deck is designed to be easily removed\* and wiped down. The headboard, footboard and side-rails are also designed with a smooth surface to facilitate cleaning and disinfection.

\*Except for the head section with X-ray cassette

Provision for Disinfection

1% sodium hypochlorite diluent can be used for disinfection purpose.

Others

Accessory Rail

Accessory rails are located on the right and left sides of the bed. These rails can be used to hang accessories such as drainage bags.



Angle Indications

Real time angle indicators for:

- Head of bed
- Trendelenburg/ Reverse trendelenburg



Other Features

- Battery Backup
- X-ray imaging deck\*

\* Applicable for selected bed models only.


Mattress (optional)



- Waterproof
- Chemical resistant (sodium hypochlorite 1.0% diluent)
- Flame retardant
- MRSA antibacterial

Accessories (optional)


- IV pole
- Oxygen cylinder holder
- Extension mattress



Reducing Shear Force

Smart Raise

Smart Raise opens the bed at the point corresponding to the sacral region to reduce shear force and prevent abdomen compression.



Retraction:

- Upward: 4-cm range
- Backward: 13-cm range

Maneuverability

Easy Grip Handle

Our easygrip handles of the head/footboards facilitate transportation. Boards are lockable to prevent accidental removal.



Address: Gate no 197, Dehu Alandi road,  
Jadhavwadi, Chikali Pune,  
Maharashtra-411062  
Mob no: +91 8377832832  
Web Site : [www.procaredmedfurn.com](http://www.procaredmedfurn.com)  
GST NO: 27ABCFP8176F1ZY

Procure

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.