

Accompanying the Recovery Journey

PLAIN GYNECOLOGICAL TABLE

SKU: PRO062

www.procaremedfurn.com

- Length: 1830 mm
- Width: 610 mm
- Height: 810 mm

KEY FEATURE


This plain gynecological examination table is made from premium CRCA steel to offer excellent built quality and strength, the table is pretreated and epoxy powder-coated for additional safety against corrosion and environmental damages and to improve the life of the product

Procure's Plain gynecological examination table comes with a 2" Premium mattress. This mattress offers excellent support and comfort to the patient during the examination and treatments.

The Single foot step stool is provided with the table which helps to easily down on the table, it helps pregnant women or elders to access the table.

This plain gynecological examination table is specially designed for gynecologist and maternity clinics, to meet their requirements to take care of women's health in their medical facility. Upgrade your hospital with Procure's Plain Gynecological Table which is the best solution to improve hospital healthcare facility.

Create a Safer Environment




Premium Quality Cushion

Powder Coated Material

CE

Procure

Address: Gate no 197, Dehu Alandi road, Jadhavwadi, Chikali Pune, Maharashtra-411062
Mob no: +91 8377832832
Web Site : www.procaremedfurn.com
GST NO: 27ABCFP8176F1ZY



Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procure