



Create a Safer Environment

PLAIN EXAMINATION TABLE

SKU: PRO070
www.procaremedfurn.com

- Length- 1830 mm
- Width- 610 mm
- Height- 530 mm

KEY FEATURE

The table is secured with a multi-layered epoxy powder coating that offers a premium look, increases durability, and resists corrosion.

This Plain OPD Table Mattress offers excellent comfort to the patient, the mattress comes with a premium leather finish and prevents water droplets and bad smells. Mattress size 2 inch

The Footstep stool helps the patient to access the table easily

This OPD table is designed for hospitals to meet their requirements and to provide a reliable and versatile solution for medical examinations and procedures. The compact design of the Plain OPD table makes it more convenient and suitable to use in different hospital wards.

Upgrade your OPD sections with our Plain Examination table specifically designed for the OPD sections, to offer excellent medical service and care to the patient.



Procare

Address: Gate no 197, Dehu Alandi road,
Jadhavwadi, Chikali Pune,
Maharashtra-411062
Mob no: +91 8377832832
Web Site : www.procaremedfurn.com
GST NO: 27ABCFP8176F1ZY

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare