



# Accompanying the Recovery Journey

## BABY CRADLE ABS

SKU: PRO023

[www.procaremedfurn.com](http://www.procaremedfurn.com)

- Width: **510 mm**
- Depth: **220 mm**
- Sleeping Basin: **580 mm**
- Length: **880 mm**
- Mattress dimensions: **740 mm L x 410 mm W x 30 mm H**
- Color Options: Available in **Blue**, and **Pink** Color options

## KEY FEATURE

This ABS baby cradle offers premium alloy for better stability and security of the cradle

The basin is built using a premium ABS material which is durable and comfortable and also enhances the safety of the baby

Its ABS body is corrosion-resistant and impact-resistant which helps to improve the lifespan of the baby cradle

A full foam mattress is provided with this ABS baby cradle to offer the best comfort to the newborn baby.

This ABS baby cradle allows you to adjust the tilt position of the basin 0° to 12°.

To maintain the smooth motion of the cradle and to help the baby to sleep easily


The cradle is equipped with anti-winding 100 mm caster wheels with brakes. these wheels are made of premium aluminum alloy which offers better grip, stability and smooth movement to the ABS Baby Cradle

Create a Safer Environment

Mattress



Bassinet



Wheel



Tilt Adjustment



# Procare

Address: Gate no 197, Dehu Alandi road,  
Jadhavwadi, Chikali Pune,  
Maharashtra-411062  
Mob no: +91 8377832832  
Web Site : [www.procaremedfurn.com](http://www.procaremedfurn.com)  
GST NO: 27ABCFP8176F1ZY

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare