



Accompanying the Recovery Journey

Create a Safer Environment

TWO FOLD STRETCHER ON TROLLEY

www.procaremedfurn.com

- Length – **2100 mm** After folding the stretcher **1050 mm**
- Width – **610 mm**
- Height – **810 mm**

KEY FEATURE


PROCARE'S TWO-SIDED FOLDABLE STRETCHER IS PRETREATED AND EPOXY POWDER-COATED WHICH OFFERS PREMIUM DURABILITY, MAKING IT RUST RESISTANT AND SHINY. THIS FOLDABLE STRETCHER IS ALSO AVAILABLE IN STAINLESS STEEL TOP AND FULL STAINLESS STEEL DESIGN. THE TOP OF THE STRETCHER IS FIXED. AND ALSO COMES WITH A NUT BOLT SETTING TO DETACH THE TOP EASILY.

THIS TWO SIDE FOLDABLE DESIGN MAKES IT EASY TO TRANSPORT SIMPLY AND REQUIRE LOW SPACE. ITS FIXED TOP AND CRCA FRAME PROVIDE STABILITY DURING DIFFERENT OPERATIONS AND TRANSPORT

PROCARE'S TWIN-FOLD STRETCHER IS EQUIPPED WITH 100 MM DIAGONAL WHEELS, ALLOWS TO MOVE STRETCHER EASILY AND KEEP IT STABLE IN TRANSPORT

THIS TROLLEY OFFERS TWO OPTIONS FOR DURABILITY, HYGIENE, AND LOOK YOU CAN CHOOSE FROM STAINLESS STEEL OR CRCA TOP OPTIONS.

PROCARE'S 2 SIDES FOLDABLE STRETCHER IS THE BEST SOLUTION FOR HOSPITALS AND CLINICS TO TRANSPORT PATIENTS ITS PREMIUM CONSTRUCTION AND DESIGN OFFER THE BEST PATIENT HANDLING AND SAFETY.



CE

Procare

Address: Gate no 197, Dehu Alandi road,
Jadhavwadi, Chikali Pune,
Maharashtra-411062
Mob no: +91 8377832832
Web Site : www.procaremedfurn.com
GST NO: 27ABCFP8176F1ZY

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare