

# Accompanying the Recovery Journey

Create a Safer Environment

## BEDSIDE SERVICE TROLLEY SS

SKU: PRO118


[www.procaremedfurn.com](http://www.procaremedfurn.com)

- Length: **460 mm**
- Width: **330 mm**
- Height: **860 mm**

### KEY FEATURE

- **MULTIPLE SHELVES:**
- PROCARE'S BEDSIDE SERVICE TROLLEY COMES WITH THREE SHELVES, THE SHELVES ARE SPACIOUS AND COME WITH THREE SIDE RAILS WHICH PREVENT ACCIDENTAL FALLS OF MEDICINES AND EQUIPMENT DURING TRANSPORT. YOU CAN ALSO BUY A TWO-SHELF VARIANT.
- **ROBUST DESIGN:**
- THIS BEDSIDE SERVICE TROLLEY IS CONSTRUCTED USING PREMIUM STAINLESS STEEL WHICH IS HIGHLY DURABLE, CORROSION-RESISTANT, AND EASY TO CLEAN. STAINLESS STEEL OFFERS A PREMIUM MIRROR FINISH. THIS TROLLEY IS ALSO AVAILABLE IN PRETREATED AND EPOXY POWDER-COATED CRCA STEEL OPTIONS.
- **SMOOTH MOBILITY:**
- THE TROLLEY COMES WITH 50 MM DIAMETER DIAGONAL WHEELS WHICH HELPS TO KEEP THE TROLLEY STEADY AND MOVE SMOOTHLY
- UPGRADE YOUR HOSPITAL WARDS WITH PROCARE'S STAINLESS STEEL BEDSIDE SERVICE TROLLEY IT IS THE PERFECT SOLUTION TO IMPROVE SERVICE DELIVERY AND PATIENT CARE.

Stainless  
Steel



Grab  
Handle




Wheel



CE

# Procare

Address: Gate no 197, Dehu Alandi road,  
Jadhavwadi, Chikali Pune,  
Maharashtra-411062  
Mob no: +91 8377832832  
Web Site : [www.procaremedfurn.com](http://www.procaremedfurn.com)  
GST NO: 27ABCFP8176F1ZY



Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare