


Removable
Stretcher



Accompanying the Recovery Journey

Create a Safer Environment

PREMIUM EMERGENCY RECOVERY TROLLEY

SKU: PR0077

www.procaremedfurn.com

Length- 2130 mm

Width- 610 mm

Height -Min. H 760 mm to Max. H 1020 mm

KEY FEATURE

THIS EMERGENCY RECOVERY TROLLEY IS CONSTRUCTED WITH PREMIUM MS PIPES MAKING IT MORE DURABLE AND STRONG.

VERSATILE SIDE RAILINGS:

FOR PATIENT SAFETY AND TO AVOID ACCIDENTAL FALLS FROM THE RECOVERY TROLLEY POWER-COATED SIDE RAILINGS ARE PROVIDED, YOU CAN CHOOSE A STAINLESS STEEL COLLAPSIBLE SIDE RAILING

SMOOTH MOBILITY:

THIS EMERGENCY RECOVERY TROLLEY COMES WITH 125 MM CASTER WHEELS, THAT PROVIDE A LOWER TURNING RADIUS TO MOVE THE TROLLEY EASILY

POSITIONS:

YOU CAN EASILY ADJUST THE BACKREST OF THE RECOVERY TROLLEY, IT ALSO PROVIDES TRENDELENBURG AND REVERSE TRENDELENBURG POSITIONS FOR BETTER COMFORT, STABILITY, AND SAFETY PRECAUTIONS DURING SURGICAL PROCEDURES.

WALL PROTECTION BUFFERS:

THE INTEGRATED SIDE BUFFERS ARE PROVIDED TO AVOID ACCIDENTAL DAMAGE, IT PROTECT FROM WALLS AND CORNERS.

Premium Quality
Cushion

Grabe
Handle

Side
Railing

Crank

Wheel



CE

Procare

Address: Gate no 197, Dehu Alandi road,
Jadhavwadi, Chikali Pune,
Maharashtra-411062

Mob no: +91 8377832832

Web Site : www.procaremedfurn.com

GST NO: 27ABCFP8176F1ZY

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare