



Accompanying the Recovery Journey

Create a Safer Environment

OT FUMIGATOR
SKU: PRO183
www.procaremedfurn.com

KEY FEATURE

- **STAINLESS STEEL CONSTRUCTION:** THE OT FUMIGATOR IS MADE FROM PREMIUM STAINLESS STEEL, ENSURING ROBUSTNESS AND RESISTANCE TO CORROSION FOR LONG-TERM USE IN HEALTHCARE FACILITIES.
- **AEROSOL PARTICLE GENERATION:** UTILIZING ADVANCED TECHNOLOGY, THE FUMIGATOR CONVERTS LIQUID FORMALIN INTO TRUE AEROSOL PARTICLES, ENSURING THOROUGH STERILIZATION OF THE OT AND ICU ENVIRONMENTS.
- **HIGH EVAPORATION RATE:** WITH AN EVAPORATION CAPACITY OF 2-2.5 LITERS PER HOUR AND A 5-LITER CAPACITY, OUR FUMIGATOR EFFICIENTLY DISPERSES STERILIZING AGENTS TO EFFECTIVELY ELIMINATE PATHOGENS AND BACTERIA.
- **WHEEL-BASED DESIGN:** MOUNTED ON A ROUNDED TRAY WITH WHEELS, THE FUMIGATOR OFFERS EASY MOBILITY AND MANEUVERABILITY, ALLOWING FOR CONVENIENT TRANSPORTATION AND PLACEMENT WITHIN MEDICAL FACILITIES.
- **VERSATILE APPLICATION:** SUITABLE FOR STERILIZING OPERATING THEATERS AND INTENSIVE CARE UNITS, OUR FUMIGATOR PROVIDES COMPREHENSIVE DISINFECTION TO MAINTAIN A HYGIENIC AND SAFE ENVIRONMENT FOR PATIENTS AND HEALTHCARE PROFESSIONALS.
- **ELECTRICAL SPECIFICATIONS:** THE FUMIGATOR OPERATES ON A 230V, 50HZ AC SUPPLY, ENSURING COMPATIBILITY WITH STANDARD ELECTRICAL SYSTEMS. ADDITIONALLY, ITS LIGHTWEIGHT DESIGN FACILITATES EFFORTLESS HANDLING AND OPERATION BY MEDICAL STAFF.

Fan



Grab
Handle



Wheel



Stainless Steel



CE

Procare

Address: Gate no 197, Dehu Alandi road,
Jadhavwadi, Chikali Pune,
Maharashtra-411062
Mob no: +91 8377832832
Web Site : www.procaremedfurn.com
GST NO: 27ABCFP8176F1ZY

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare