


Procare



Accompanying the Recovery Journey

Create a Safer Environment

PREMIUM CARDIAC TABLE

SKU: PRO112

www.procaremedfurn.com

- Length: **410 mm**
- Width: **560 mm**
- Height Adjustment: **890 mm to 1400 mm**

KEY FEATURE

CARDIAC TABLE USES

YOU CAN USE A CARDIAC TABLE TO PERFORM THE FOLLOWING ACTIVITIES

- 1.THERAPEUTIC USE IN CARDIAC SURGERY
- 2.IT CAN BE USED TO KEEP PRESCRIPTION DRUGS AND MEDICAL DEVICES.
- 3.EATING A FOOD
- 4.WRITING AND READING
- 5.TO USE LAPTOP OR TABLETS
- 6.MONITORING PATIENT
- 7.ASSISTING THERAPEUTIC ACTIVITIES. ETC.

CARDIAC TABLE: SPECIFICATIONS AND FEATURES

HEIGHT ADJUSTMENT:

THIS PREMIUM CARDIAC TABLE COMES WITH A HEIGHT ADJUSTMENT FUNCTION, IN WHICH YOU CAN EASILY ADJUST THE HEIGHT OF THE CARDIAC TABLE AS YOU REQUIRE. THE HEIGHT ADJUSTMENT RANGES FROM 890 MM TO 1400 MM

PREMIUM ABS TOP:

PROCARE'S CARDIAC TABLE OFFERS DIFFERENT TOP OPTIONS SUCH AS ABS TOPS, WOODENTOPS, AND STAINLESS STEEL TOPS TO COMPLETE THE DIFFERENT REQUIREMENTS OF THE HOSPITAL IN DIFFERENT TREATMENTS, EXAMINATIONS, AND DIAGNOSTIC PROCEDURES, AND TO PROVIDE THE BEST PLATFORM FOR CARDIAC ACTIVITIES.

PREMIUM CONSTRUCTION:

THE CARDIAC TABLE IS CONSTRUCTED FROM HIGH-QUALITY ALUMINUM ALLOY. THE BASE OF THE TABLE IS MADE FROM PREMIUM CRCA STEEL AND COATED WITH EPOXY POWDER WHICH OFFERS THE BEST DURABILITY AND STRENGTH TO THE TABLE AND MAKES IT CORROSION-RESISTANT. IT HELPS TO USE A TABLE IN DIFFERENT CONDITIONS.

PREMIUM WHEELS:

Height
Adjustable

ABS Top

Wheel



Procare

CE


Procare

Address: Gate no 197, Dehu Alandi road,
Jadhavwadi, Chikali Pune,
Maharashtra-411062

Mob no: +91 8377832832

Web Site : www.procaremedfurn.com

GST NO: 27ABCFP8176F1ZY



Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare