

# Accompanying the Recovery Journey

Create a Safer Environment

PROCARE RADIANT  
WARMER STAND-  
ALONE  
[www.procaremedfurn.com](http://www.procaremedfurn.com)

## KEY FEATURE

- Micro Processor Based servo Controller for precise temperature regulation.
- Multiple Modes of Operation including Servo skin, AIR, Pre-Warm, and Manual, catering to diverse clinical needs.
- Two digital displays enable accurate monitoring of baby's temperature and set temperature.
- User-friendly soft touch buttons facilitate intuitive operation.
- Convenient °F and °C Converter for flexibility in temperature units.
- The Baby Feed Reminder feature ensures timely feeding schedules are maintained.
- Integrated Stand for Monitor and I.V. provides easy access to essential equipment.
- Swiveling top allows for X-ray procedures without disrupting the infant's position.
- LED Illumination ensures clear observation of the baby's condition.
- Safety alarm with high-temperature cut-off offers additional protection for newborns.

IV Stand

Adjustable  
Knob


Wheel



CE

**Procare**

Address: Gate no 197, Dehu Alandi road,  
Jadhavwadi, Chikali Pune,  
Maharashtra-411062  
Mob no: +91 8377832832  
Web Site : [www.procaremedfurn.com](http://www.procaremedfurn.com)  
GST NO: 27ABCFP8176F1ZY



Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

**Procare**