

# Accompanying the Recovery Journey

SUCTION MACHINE WITH POLYCARBONATE JAR  
SKU: PRO191  
[www.procaredmedfurn.com](http://www.procaredmedfurn.com)

Create a Safer Environment

## KEY FEATURE

- BETTER CAPACITY:**  
THIS MACHINE COMES WITH A STRONG 700 ± 10 MM HG SUCTION CAPACITY AND IT IS ADJUSTABLE TO 20-30 LMP, TO ENSURE A EFFICIENT SUCTION.
- VACUUM MONITORING:**  
YOU CAN MONITOR AND CONTROL SUCTION LEVELS WITH THE HELP OF A 2-INCH BCPM VACUUM GAUGE. THE GAUGE IS CALIBRATED BETWEEN 0-760 MM HG.
- POLYCARBONATE JARS**  
THE SUCTION MACHINE IS EQUIPPED WITH A SET OF TWO POLYCARBONATE JARS WITH SILICON SEALING AND BLUE COLORED ABS LIDS, SPECIFICALLY DESIGNED TO ENSURE SAFE AND HYGIENIC FLUID CONTAINMENT.
- STRONG PUMP:**  
THE MACHINE OFFERS A 1/4 HORSEPOWER ELECTRIC MOTOR WITH AN OIL-IMMERSED PISTON, TO PROVIDE POWERFUL AND CONSTANT SUCTIONS
- USER-FRIENDLY DESIGN:**  
THE MACHINE COMES WITH A STRONG CABINET MADE FROM POWDER-COATED M.S. THE MACHINE IS MOUNTED ON RUBBER BUSHINGS FOR STABILITY AND HAS A HANDLE FOR BETTER HOLDING.
- EASY TO TRANSPORT:**  
THIS SUCTION MACHINE IS COMPACT AND EASY TO TRANSPORT THE WEIGHT OF THE MACHINE IS 19 KG AND THE DIMENSIONS ARE 19" X 12" X 29".



Jar



Grab Handle




Wheel



# Procure

Address: Gate no 197, Dehu Alandi road, Jadhavwadi, Chikali Pune, Maharashtra-411062  
Mob no: +91 8377832832  
Web Site : [www.procaredmedfurn.com](http://www.procaredmedfurn.com)  
GST NO: 27ABCFP8176F1ZY



Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

