

Accompanying the Recovery Journey

Create a Safer Environment

AMBULANCE STRETCHER

SKU: PR0080

www.procaremedfurn.com

- Normal Dimensions: Length-**190 cm**, Width-**55 cm**, Height-**90 cm**
- Lowest Height Position of the Ambulance Stretcher: Length-**190**, Width-**55**, Height- **25 cm**
- Maximum Height Position of the Foldable Ambulance Stretcher: Length- **190 cm**, Width-**55cm**, Height- **86 cm**
- Loading Capacity- 150 kg
- Construction- aluminum
- Wheel diameter: 125 mm

KEY FEATURE

- **COMPLETE FOLDABLE**
- **SIDE RAILS**
- **HOLDING BELTS**
- **BACKREST**
- **NOTE: THE ACTUAL PRODUCT MAY BE DIFFERENT FROM THE IMAGE; IT IS ONLY MEANT TO BE USED AS A GUIDE.**

THE HEIGHT OF THIS AMBULANCE STRETCHER IS ADJUSTABLE TO FIT IN DIFFERENT EMERGENCY CONDITIONS. THE HEIGHT OF THIS AMBULANCE STRETCHER IS RANGING FROM 25 CM TO 86 CM.

THIS FOLDABLE AMBULANCE STRETCHER OFFERS 150 KG OF LOADING CAPACITY WITH A PREMIUM ALUMINUM CONSTRUCTION TO PROVIDE BEST PERFORMANCE AND COMFORT TO THE PATIENT

THE SIDE RAILING HELPS TO TRANSPORT THE PATIENT SAFELY TO THE AMBULANCE BY PREVENTING ACCIDENTAL FALLS.

THE HIGHLY STRENGTH BELTS HELP TO HOLD THE PATIENT ON THE TROLLEY, ONE BELT IS PROVIDED FOR THE CHEST AREA, AND ANOTHER FOR THE LEGS TO PREVENT MOVEMENT OF THE PATIENT DURING TRANSPORT

Safety
Belt




Railing



Wheel




CE



Procare

Address: Gate no 197, Dehu Alandi road,
Jadhavwadi, Chikali Pune,
Maharashtra-411062
Mob no: +91 8377832832
Web Site : www.procaremedfurn.com
GST NO: 27ABCFP8176F1ZY



Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare