

Accompanying the Recovery Journey

CHIROPRACTIC DROP FIVE TABLE

www.procaremedfurn.com

Minimum Height of the table – 21

Maximum Height: 32 inches is the maximum height.

Measurement in length: 74 inches

Breadth: 26 inches

Max. Load: 180 kg

Weight of the Table: 72 kg net weight

Warranty of the Product: One-year warranty

KEY FEATURE

Our Chiropractic tables come with five manual drops for therapy including drops at the Cervical, Thoracic, Lumbar, and Pelvic Regions.

Height Adjustable:

You can easily adjust the height of the chiropractic drop five tables the height ranges from 21 to 32 inches for better positioning and easy operations.

Robust Construction:

The table is constructed on premium CRCA steel and the body is powder coated for extra safety and finishing, which makes it rust-resistant and easy to lean.

High Load Capacity:


Our chiropractic drop-five table has a massive weight capacity the maximum weight capacity of the table is 180 kg.

Compact Design:

Its unique design makes it more adaptable and easy to adjust position and move between the wards with its retractable wheels.

Ideal Dimensions:

This table has 74 Inch Length and 26-inch Width, which offers enough room for the patients to feel comfortable




Create a Safer Environment

CE

Procare

Address: Gate no 197, Dehu Alandi road,
Jadhavwadi, Chikali Pune,
Maharashtra-411062
Mob no: +91 8377832832
Web Site : www.procaremedfurn.com
GST NO: 27ABCFP8176F1ZY



Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare