

# Accompanying the Recovery Journey

Create a Safer Environment

## DELUXE ATTENDANT BED

**SKU: PR0017**

**[www.procaremedfurn.com](http://www.procaremedfurn.com)**

Dimensions of the attendant beds are.

- Length: **1830 mm**
- Width: **620 mm**
- Height: **460 mm**

## KEY FEATURE

### ROBUST DESIGN:

THIS DELUXE ATTENDANT BED IS MADE FROM PREMIUM CRCA STEEL, ITS PREMIUM STEEL OFFERS TOUGHNESS AND DEPENDABILITY

### SUPERIOR SURFACE:


THE ATTENDANT BED WAS PRETREATED AND EPOXY POWDER-COATED TO CREATE A SMOOTH, CORROSION-RESISTANT, AND LONG-LASTING SURFACE

### STRONG SUPPORT:

THIS DELUXE ATTENDANT COMES WITH METAL PLANS THAT HAVE A STRONG WEIGHT CAPACITY, TO GIVE SUPPORT TO THE PATIENT




CE



# Procare

Address: Gate no 197, Dehu Alandi road,  
Jadhavwadi, Chikali Pune,  
Maharashtra-411062  
Mob no: +91 8377832832  
Web Site : [www.procaremedfurn.com](http://www.procaremedfurn.com)  
GST NO: 27ABCFP8176F1ZY



Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

**Procare**