

# Accompanying the Recovery Journey



Storage space For Medicine.

## MANUAL EMERGENCY RECOVERY TROLLEY

[www.procaremedfurn.com](http://www.procaremedfurn.com)

- The Length of the Manual Recovery Trolley- **2100 mm**,
- The width of the Recovery Trolley –**710 mm**,
- Min. Height – **650 mm**,
- Max. Height- **950 mm**
- Pretreated and Epoxy powder-coated frame.
- Different positions adjustable manually using screw- Height, Backrest, Trendelenburg, and reverse Trendelenburg positions

### KEY FEATURE

- OXYGEN CYLINDER CAGE
- I.V. ROD
- SIDE RAILINGS WITH SWING-AWAY DESIGN (PC)
- 125MM-DIAMETER, HEAVY-DUTY CASTOR WHEELS WITH BRAKES
- SIDE BUFFERS AT CORNERS
- 

PROCARE'S MANUAL EMERGENCY RECOVERY TROLLEY IS A GOOD COMBINATION OF FEATURES AND COMFORTS ITS PREMIUM SCREW MECHANISM HELPS TO ADJUST THE HEIGHT AND OTHER POSITIONS ACCURATELY.


THE FRAME OF THIS MANUAL EMERGENCY RECOVERY TROLLEY IS PRETREATED AND EPOXY POWDER-COATED TO RESIST CORROSION AND PROVIDE LONG LIFE AND DURABILITY. THIS MAKES THE TROLLEY PERFECT TO USE IN HARSH CONDITIONS

. THE POWDER-COATED SWING-AWAY SIDE RAILING PROVIDES MORE SAFETY AGAINST UNINTENTIONAL FALLS OF THE PATIENT DURING TRANSPORT

IT IS ANOTHER ESSENTIAL PART OF ANY EMERGENCY INSTRUMENT, THIS OXYGEN CYLINDER CAGE HELPS TO CARRY A CYLINDER WITH A TROLLEY TO PROVIDE LIFE SUPPORT TO THE PATIENT.



X-ray permeable Sheet



Cage for Oxygen Cylinder provision



Stainless Steel Swing-Away Side Railings for safety



Grab Hndle

Four Heavy-Duty Swivel Castors (125mm Diameter) with Brakes

X-Ray Cassette Tray Under Stretcher Top



Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

# Procare

Address: Gate no 197, Dehu Alandi road,  
Jadhavwadi, Chikali Pune,  
Maharashtra-411062  
Mob no: +91 8377832832  
Web Site : [www.procaremedfurn.com](http://www.procaremedfurn.com)  
GST NO: 27ABCFP8176F1ZY

Procare