

Create a Safer Environment

## Accompanying the Recovery Journey

### EXAMINATION COUCH 3 DRAWER 3 CABINETS [www.proccaremedfurn.com](http://www.proccaremedfurn.com)

- Length: 1830 mm
- Width: 610 mm
- Height: 810 mm
- Backrest: 0° - 70°

#### KEY FEATURE

This Examination couch offers an adjustable backrest for extra comfort you can adjust it up to 70°.

This examination couch with a backrest is built from CRCA steel, which provides more durability and longer life. The Couch is pretreated and Epoxy Powder Coated, which resists corrosion, and improves product durability to regular use.

The couch with a backrest comes with a 2-inch mattress offers excellent comfort and support to the patient during examination procedures.

The Examination couch with backrest comes with integrated three drawers, three cabinets, and one BP Apparatus tray. for a spacious storage solution to store emergency medicines.

A Single footstep stool is provided with the Examination Couch.

Adjustable Backrest →

Premium Quality Cushion ↙




Proccare

CE

**Proccare**

Address: Gate no 197, Dehu Alandi road,  
Jadhavwadi, Chikali Pune,  
Maharashtra-411062  
Mob no: +91 8377832832  
Web Site : [www.proccaremedfurn.com](http://www.proccaremedfurn.com)  
GST NO: 27ABCFP8176F1ZY



Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Proccare