

## Mattress



# Accompanying the Recovery Journey

## Create a Safer Environment

[illegible]

## KEY FEATURE

To complete different medical needs and obtain a perfect position you can easily raise or lower the bed's height from (Approx.) 385mm to 710mm


To improve patient for patients who need knee rests, by offering a comfortable nee rest angle of 30° to 35°.

The backrest is an essential function for the patient's comfort, This bed offers up to 75° backrest angle for better support and comfort.

This Manual ICU Bed with Railings provides a 12° - 14° angle Trendelenburg position.

The Reverse Trendelenburg angle of this bed is 8°-10° which helps in different medical treatments.

You can add These extra accessories to your Manual five-function ICU Bed With Collapsible Railings to make your bed more adaptable.



Collapsible  
Railing

# Procure

**Address: Gate no 197, Dehu Alandi road,  
Jadhavwadi, Chikali Pune,  
Maharashtra-411062**  
**Mob no: +91 8377832832**  
**Web Site : [www.procuremedfurn.com](http://www.procuremedfurn.com)**  
**GST NO: 27ABCFP8176F1ZY**

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

ProCare