



Create a Safer Environment

# Accompanying the Recovery Journey

## FOWLER BED MANUAL WITH ABS HEAD AND FOOTBOARD

[www.procaremedfurn.com](http://www.procaremedfurn.com)

- length of the Manual Fowler Bed: **2240 mm**
- Width of the Manual Fowler Bed: **890 mm**
- Height of the Manual Fowler Bed: **560 mm** in height (without mattress)

### KEY FEATURE

- This Fowler bed with ABS Head and Footboard is manually adjustable, the backrest of the bed is adjustable up to 75° and the knee rest of the bed is adjustable between 30° – 35° easily using the screw mechanism.
- Screw mechanism Cranks are provided to control positions manually.

This manual Fowler bed is constructed on Premium Powder coated CRCA steel, Which is strong, and rust-resistant.

A high-quality metal sheet is used to construct the bed which supports up to 125 kg weight.

These premium quality bows help to look bed premium and also increase the safety of the bed

Easily visible and accessible I.V. rod provision.



crank

ABS Panel

Raichet  
System

CE

# Procare

Address: Gate no 197, Dehu Alandi road,  
Jadhavwadi, Chikali Pune,  
Maharashtra-411062  
Mob no: +91 8377832832  
Web Site : [www.procaremedfurn.com](http://www.procaremedfurn.com)  
GST NO: 27ABCFP8176F1ZY

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare