

Accompanying the Recovery Journey

Create a Safer Environment

FULL GYNECOLOGICAL EXAMINATION COUCH

SKU: PRO064

www.procaremedfurn.com

- Length: 1830 mm
- Width: 610 mm
- Height: 810 mm

KEY FEATURE

This examination couch is made from premium CRCA steel material. this examination couch offers stability and a long life. This couch is layered with an epoxy powder coating to protect against rust and scratches.

This examination couch offers you a spacious storage solution, three drawers, four cabinets, and One BP Apparatus tray, doctors can store all medicines or equipment required in emergencies.

A smooth and comfortable 2-inch mattress is provided with the couch to improve the comforts of the couch, this mattress is easy to clean and water droplet resistant, it also resists bad smells and helps to maintain hygiene.

A single-footstep stool helps the patient to in and out of the couch easily.

This full gynecological examination couch is available in different colors Ex. Blue, Green, grey etc.



Lithotomy Rod

Steap Stool


BP Tray

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare

Procare

Address: Gate no 197, Dehu Alandi road, Jadhavwadi, Chikali Pune, Maharashtra-411062
Mob no: +91 8377832832
Web Site : www.procaremedfurn.com
GST NO: 27ABCFP8176F1ZY



CE

