



Accompanying the Recovery Journey

ONE DRAWER THREE CABINET GYNECOLOGICAL EXAMINATION COUCH

SKU: PRO066

www.procaremedfurn.com

- Length: 1830 mm
- Width: 610 mm
- Height: 810 mm

KEY FEATURE

The whole body of this gynecological couch is made from CRCA steel and coated with an epoxy powder. Its coating makes this couch corrosion-resistant and scratch-resistant and improves the durability of the product.

To store emergency medicine and equipment multiple storage options are provided with the couch, One drawer, Three cabinets, and one BP apparatus tray are provided

To provide extra comfort to the patient 2" mattress is provided with the couch. this mattress comes with a premium leather finish

Create a Safer Environment



Procare

Address: Gate no 197, Dehu Alandi road,
Jadhavwadi, Chikali Pune,
Maharashtra-411062
Mob no: +91 8377832832
Web Site : www.procaremedfurn.com
GST NO: 27ABCFP8176F1ZY

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare