



HYDRAULIC LABOUR TABLE

SKU: PRO068
www.procaremedfurn.com

- Length: 1880 mm
- Width: 900 mm
- Height: 600 mm to 800 mm

Accompanying the Recovery Journey

Create a Safer Environment

KEY FEATURE

Procure's Hydraulic Labour Table comes with a foot-operated Hydraulic mechanism, this mechanism helps medical professionals to adjust the height of the table easily and also helps to comfort during childbirth.

Procure's hydraulic Labour table offers a backrest, Trendelenburg, and Revere Trendelenburg position to offer the best comfort and positioning. The top of the Labour table offers a stainless steel three-section top which is easy to clean

This Hydraulically operated trolley is made from premium Mild Steel that offers excellent durability and long life.

This Hydraulically operated trolley is made from premium Mild Steel that offers excellent durability and long life.

The hydraulic Labour table is equipped with 100 mm caster wheels, these wheels help to move the trolley easily in emergencies.

IV Stand



Leg Support



Wheel



Premium Quality Cushion



Railing



Procure

Address: Gate no 197, Dehu Alandi road,
Jadhavwadi, Chikali Pune,
Maharashtra-411062
Mob no: +91 8377832832
Web Site : www.procaremedfurn.com
GST NO: 27ABCFP8176F1ZY

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procure