

# Accompanying the Recovery Journey

## LAB SERIES PADDED BLOOD DRAWING CHAIR

[www.procaremedfurn.com](http://www.procaremedfurn.com)

- Length: **610 mm**
- Width: **610 mm**
- Height: **960 mm**

### KEY FEATURE

THIS LAB SERIES PADDED BLOOD COLLECTION CHAIR IS CONSTRUCTED USING CRCA STEEL AND LAYERED WITH EPOXY POWDER TO OFFER THE BEST DURABILITY, A SHINY LOOK, AND A RUST-RESISTANT SURFACE.

#### PADDED CUSHIONED ARMREST:

THE HEIGHT ADJUSTABLE PADDED ARMREST PROVIDED WITH THE BLOOD-DRAWING CHAIR MAKES BLOOD-DRAWING OPERATIONS MORE COMFORTABLE THE PADDED ARMREST OFFERS SUPPORT TO THE ARMS BY PROVIDING ACCURATE ADJUSTMENTS AND POSITIONS

#### LOADING CAPACITY:

THIS LAB SERIES PADDED BLOOD DRAWING CHAIR IS ADAPTABLE TO DIFFERENT HOSPITALS AND CLINICS IT COMES WITH A 190 KG WEIGHT CAPACITY.

#### PREMIUM CUSHION:

THE ANTIMICROBIAL PREMIUM PU LEATHER FINISHED CUSHION OFFERS HYGIENE, COMFORT, AND PREMIUM LOOK TO THE CHAIR.


Create a Safer Environment



Arm Rest  
Up & Down




Adjustable  
Hand Rest



Premium Quality  
Cushion




CE



# Procare

Address: Gate no 197, Dehu Alandi road,  
Jadhavwadi, Chikali Pune,  
Maharashtra-411062  
Mob no: +91 8377832832  
Web Site : [www.procaremedfurn.com](http://www.procaremedfurn.com)  
GST NO: 27ABCFP8176F1ZY



Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare