



Accompanying the Recovery Journey

Create a Safer Environment

BABY CRADLE POWDER COATED

SKU: PRO024

www.procaremedfurn.com

- Length: 920 mm
- Width: 530 mm
- Height: 960 mm

KEY FEATURE

This cradle is Pretreated and epoxy powder coated, which makes it shiny, easy to clean, and corrosion-resistant and helps to maintain the hygiene of the baby and cradle.

You can easily detach the cradle from the frame if you wish to carry a cradle without the frame

This PC Baby cradle offers a silent and smooth swing which helps the baby fall asleep without any disturbance

The depth of the cradle prevents the accidental fall of the baby and ensures more safety.

Powder Coated



CE

Procare

Address: Gate no 197, Dehu Alandi road,
Jadhavwadi, Chikali Pune,
Maharashtra-411062

Mob no: +91 8377832832

Web Site : www.procaremedfurn.com

GST NO: 27ABCFP8176F1ZY

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare