



SEMI-FOWLER PEDIATRIC BED WITH ABS PANEL

www.procaremedfurn.com

- Length: 1575 mm
- Width: 860 mm
- Height: 570 mm
- Backrest – 70°

KEY FEATURE

introducing Procare's best semi fowler pediatric bed designed and built to provide maximum comfort and safety to the kids while their medical treatments


This pediatric bed offers a backrest for the support and comfort of the kids. This semi fowler bed offers a 70° backrest

The kid's bed is pre-treated and epoxy powder-coated, which provides a corrosion-resistant surface and extends the durability of the bed, making it smooth and shiny.

This semi-fowler bed is made from complete CRCA steel which increases the strength and durability of the semi-fowler pediatric bed and improves the safety of the kids.


The pediatric bed is equipped with full-length drop-side railings to improve the kid's safety and prevent accidental falls

Two premium stainless steel I.V. rods are provided with the bed for easy access and complete kids' care.



Accompanying the Recovery Journey

Create a Safer Environment



CE

Procare

Address: Gate no 197, Dehu Alandi road,
Jadhavwadi, Chikali Pune,
Maharashtra-411062
Mob no: +91 8377832832
Web Site : www.procaremedfurn.com
GST NO: 27ABCFP8176F1ZY

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare