

Procare

Accompanying the Recovery Journey

Create a Safer Environment

PLAIN BACKREST EXAMINATION TABLE

www.procaremedfurn.com

- Length: 1830 mm
- Height: 810 mm
- Width: 610 mm

KEY FEATURE

The examination table is made from pretreated and epoxy powder-coated premium CRCA steel, which offers the best durability, and corrosion-resistant smooth surface.

This Plain OPD Table comes with a backrest adjustment, you can adjust the backrest up to 75° to offer the best comfort during the examination.

A 2-inch premium Mattress is provided with the Plain OPD examination table with a backrest, this mattress comes with a premium leather finish, which prevents bad smells, and water droplets and offers a clean surface to maintain the hygiene of the table

A Single foot step stool helps the patient to use the table easily.

Adjustable
Backrest

Premium Quality
Cushion

Raichet
System

Procare

Address: Gate no 197, Dehu Alandi road,
Jadhavwadi, Chikali Pune,
Maharashtra-411062
Mob no: +91 8377832832
Web Site : www.procaremedfurn.com
GST NO: 27ABCFP8176F1ZY

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare