





Accompanying the Recovery Journey

Create a Safer Environment

**C-ARM HYDRAULIC
OT TABLE**
SKU: PROOT002
www.procaremedfurn.com

- **Top Dimensions of the Table** – Length 1980mm X Width 530mm
- **Height Adjustment** -730mm to 1030mm
- **Reverse Trendelenburg** – 20° / 25°
- **Lateral Tilt** – 20° / 20°
- **Kidney Elevator** – 130mm
- **Back Rest (UP and Down)** – 80° / 35°
- **Leg Rest (Up and Down)**- 10° / 90°
- **Head Rest (Up and Down)**- 45° / 90°
- **Patient Weight Capacity** – 200 Kgs

KEY FEATURE

effortless positioning and adjusting functions of the OT Table can be operated Mechanically and Hydraulically.

designed with slanted pillars to provide for clear C-Arm moments, improving image performance.

The Table has come with hygienic and long-lasting stainless steel fittings that are simple to maintain.

The pelvic region, head, and legs can all be detached enabling simple configuration and flexible patient accommodations.

The bed's unique flexible design allows you to change the head and leg sections to provide better support to the patient.

The inbuilt kidney bridge of this OT table helps in different medical treatments and urological procedures.

Adjustable
Head Rest

Shoulder
Support

Waist
Support

Height
Adjustable

Premium Quality
Cushion

Leg Support

CE

Procare

Address: Gate no 197, Dehu Alandi road,
Jadhavwadi, Chikali Pune,
Maharashtra-411062
Mob no: +91 8377832832
Web Site : www.procaremedfurn.com
GST NO: 27ABCFP8176F1ZY

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare