

Accompanying the Recovery Journey

CHIROPRACTIC DROP FOUR TABLE

SKU: PROT035

www.procaremedfurn.com

You can take advantage of its large surface area, the surface area of this drop-four table is 21 Inches which offers more space to patient and treat patients. The height of the table is manually adjustable from 18 to 25 inches (46 to 63.5 cm). It makes it more usable for different treatment methods.

KEY FEATURE

Its high-density foam provides more comfort and support to the patients in different treatments.

Advanced Headpiece:

Procare's Chiropractic Drop Four Table has a multidirectional headpiece with a split cushion and armrest, and they can tilt and elevate with an option of exact adjustment.

Elevating parts:

This drop 4 table offers more functions for a precise treatment it's designed with a leg extension which helps in knee break, lumbar, and pelvic parts

Foot-Activated Drops:

The pedals are provided to the table to control and active the drop section quickly and easily

Robust Design:

This chiropractic drop-four table is constructed using premium steel, which offers long life, and durability to the product

Specialized Drop Sections:

Procare's Drop Four table offers different drop sections with elevating and drop-only options for the cervical, upper thoracic, thoracic, and lumbar regions

Customizable:

Procare can make customization as you require in the product.


Create a Safer Environment



CE

Procare

Address: Gate no 197, Dehu Alandi road,
Jadhavwadi, Chikali Pune,
Maharashtra-411062
Mob no: +91 8377832832
Web Site : www.procaremedfurn.com
GST NO: 27ABCFP8176F1ZY



Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare