

Accompanying the Recovery Journey

Create a Safer Environment

KEY FEATURE

For safety and to keep the patient steady tilt table has three broad straps, Straps help to hold the patient during procedures.

The tables offer manual functionality its worm gear system provides accurate control to raise patients into a standing posture.

The tilt Table is designed to complete all patient needs and treatment requirements, This table is calibrated to handle a tilt range of 0 to 85 kg

The table provides better space for patients with its padded surface. the table is 198 cm Long, and 61 cm Wide and 80 cm in Height.


This Tilt Table is supported by 35 mm square steel tube legs that provide more stability and safety. It has a Powder Coated with a premium finishing it makes table Rust resistant and Durable.

The tilt Table can be moved and locked with its four lockable 10 cm diameter casters. You can easily move table using these casters.

CHIROPRACTIC MANUAL TILT TABLE

SKU: PROT030

www.procaremedfurn.com



Procare

Address: Gate no 197, Dehu Alandi road,
Jadhavwadi, Chikali Pune,
Maharashtra-411062
Mob no: +91 8377832832
Web Site : www.procaremedfurn.com
GST NO: 27ABCFP8176F1ZY

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare