





Accompanying the Recovery Journey

DERMATOLOGY CHAIR

SKU: PROOT008

www.procaremedfurn.com

• Common Use:	Hospital Bed	
• Used Material:	Sponge, PVC/PU, Synthetic Leather	
• Dimensions of the Bed:	The dimensions are L 78" x W 22" x H 26"-36". The loading capacity is 200 kg.	
• Remote Control:	Provided	
• Motors:	Three motors	
• Pillow:	Individual	
• Options for Color:	Diverse	
• Uses:	Derma Clinic, Hair Transplant, Skin Cosmetology, Hospital	
• Features:	Battery Backup – Optional=	
	—	

KEY FEATURE

Introducing our best-quality Dermatology chair made to provide the best comfort to the patient and doctor. This dermatology chair is equipped with high-quality synthetic leather (PVC/PU) and sponge to make it more comfortable, soft, and durable for regular use. The length of the derma chair is 78", Width is 22", and 36" in Height. This derma chair comes with a 200 Kg capacity with height adjustment feature.

This chair is suitable for different procedures like skin cosmetology, derma clinics, and hair transplant procedures. To provide simple and easy adjustment this chair provides three motors a remote control and a battery backup option. This chair is available in various color options.

Create a Safer Environment



Adjustable Backrest


Adjustable Knee Rest

Height Adjustable

CE

Procare

Address: Gate no 197, Dehu Alandi road, Jadhavwadi, Chikali Pune, Maharashtra-411062
Mob no: +91 8377832832
Web Site : www.procaremedfurn.com
GST NO: 27ABCFP8176F1ZY



Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare