

Accompanying the Recovery Journey

Create a Safer Environment

CHIROPRACTIC ELECTRIC TILT TABLE

SKU: PROT031

www.procaremedfurn.com

height from 20 to 36 inches,
support up to 200 kg
24 by 72 Inches,

KEY FEATURE

This table allows chiropractors to adjust tables as per their requirements and procedures its premium electric motor smoothly raises or lowers the height from 20 to 36 inches,


The table is remote-operated, you can control the tilting range up to Ninety degrees (90°) very smoothly.

You can move the electric tilt table anywhere with the help of its lockable Casters. It increases the uses ad flexibility of the electric tilt table

To provide more comfort to the patient while treatment the table has a 3- inch foam cushion, Cushion is available in different colors.

This Chiropractic Electric Tilt table is built using a premium heavy duty steel frame which support up to 200 kg


At24 by 72 Inches, the table offers best room for patient of all shapes and sizes. and provide comfort during therapy.



CE

Procare

Address: Gate no 197, Dehu Alandi road,
Jadhavwadi, Chikali Pune,
Maharashtra-411062
Mob no: +91 8377832832
Web Site : www.procaremedfurn.com
GST NO: 27ABCFP8176F1ZY



Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare