





# ELECTRIC LABOUR BED

Elevate Comfort with Our Premium Hospital Bed



## Accompanying the Recovery Journey

PROCARE'S ELECTRIC DELIVERY ROOM LABOUR TABLE OFFERS DIFFERENT POSITIONS TO PROVIDE EXCELLENT COMFORT AND ADJUSTMENTS. THE BED OFFERS BACKREST, TRENDLENBURG, REVERSE TRENDLENBURG, AND HEIGHT ADJUSTMENT POSITIONS



Create a Safer Environment

[www.procaredmedfurn.com](http://www.procaredmedfurn.com)

	PRO027 (ELECTRIC LABOUR BED)	
No. of Motor	0 motor	
Total Length	1830 mm	
Total Width*1	610 mm	
Total Height	385 mm to 710 mm	
Backrest	70° degree	
Trendelenburg	0°-12°	
Mattress Deck	130 kg	
Safe Working Load*3		

CE


\*1 When the side-rail is raised.  
\*2 Indicates the height from floor to the mattress deck surface.  
\*3 Maximum load with which the bed can be operated. (The total load of user's weight and incidental items including the mattress and options)

**Procare**

Address: Gate no 197, Dehu Alandi road,  
Jadhavwadi, Chikali Pune,  
Maharashtra-411062  
Mob no: +91 8377832832  
Web Site : [www.procaredmedfurn.com](http://www.procaredmedfurn.com)  
GST NO: 27ABCFP8176F1ZY

## KEY FEATURE

The ABS side railings are provided for extra safety, The bed is constructed using premium Mild steel which offers extra durability of the bed



Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare