



# Accompanying the Recovery Journey

Create a Safer Environment

## ELECTRIC OPERATING TABLE

SKU: PROOT005

[www.procaremedfurn.com](http://www.procaremedfurn.com)

- **Top Dimension**- Length 2150 mm X Width 535 mm
- **Height Adjustment**- 750 mm to 1050 mm
- **Trendelenburg ( Reverse )** – 20° / 25°
- **Lateral Tilt**- 20° / 20°
- **Kidney Elevator**- 130 mm
- **Back Rest ( Up and Down )**- 80° / 15°
- **Leg Rest ( Up and Down )**- 10° / 90°
- **Head Rest ( UP and Down )**- 45° / 90°
- **Mains Voltage Capacity**- AC 220V±10%50Hz
- **Patient Weight Capacity**- 200kgs

### KEY FEATURE

Introducing our best surgical technology- the Electric Operating Table. Specially made to fulfill all requirements of the modern challenges in operation theatre.

This electric OT Table is completely electric which improves the handling of the OT Table. Which guarantees simple and easy adjustment for a range of surgical positions.

Our table comes with five five-function radio translucent table tops which allows for throughout imagining without requiring positioning, which improves surgical accuracy and productivity.

This table's strong foundation and eccentric pillar construction offer strong stability and support even during a hard procedures.

The head, and leg sections of the table are manually exchangeable, You can easily change these sections of the table

This Electric OT Table comes with a Cable remote control, This remote has all functions displayed with a button. every position of the table can operate using the remote,

Adjustable  
Head Rest

Shoulder  
Support

Waist  
Support

Premium Quality  
Cushion

Height  
Adjustable

Leg Support



CE

# Procare

Address: Gate no 197, Dehu Alandi road,  
Jadhavwadi, Chikali Pune,  
Maharashtra-411062  
Mob no: +91 8377832832  
Web Site : [www.procaremedfurn.com](http://www.procaremedfurn.com)  
GST NO: 27ABCF8176F1ZY

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare