

Create a Safer Environment

## Accompanying the Recovery Journey

### EXAMINATION COUCH 2 DRAWER 2 CABINETS

**SKU: PRO074**

**[www.procaremedfurn.com](http://www.procaremedfurn.com)**

- Length: 1830 mm
- Width: 610 mm
- Height: 810 mm

### KEY FEATURE

This examination couch is the perfect OPD solution for clinics and hospitals to provide excellent comfort and convenience.

This compact examination couch is made from CRCA steel to improve product strength and durability and make this product corrosion-resistant.

This examination couch offers the best storage options for storing emergency medicines and equipment it comes with Two drawers, Two cabinets, and One BP apparatus. It helps keep your things organized and easily accessible.

This couch comes with a multilayer epoxy powder coating which offers a smooth and shiny surface, that resists corrosion and keeps the product safe from water and dust.


The 2-inch premium mattress is provided with a couch to offer the best support and comfort to the patient during an examination. The mattress resists bad smell, water droplets, and sweat

The footstep accessory is provided with the examination couch it helps the patient to get onto the examination couch easily. It is an essential accessory for elders, kids, and pregnant women.



**Procare**

Address: Gate no 197, Dehu Alandi road,  
Jadhavwadi, Chikali Pune,  
Maharashtra-411062  
Mob no: +91 8377832832  
Web Site : [www.procaremedfurn.com](http://www.procaremedfurn.com)  
GST NO: 27ABCFP8176F1ZY



Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

**Procare**