

# Accompanying the Recovery Journey

Create a Safer Environment

## FOOT SUCTION MACHINE M/S

SKU: PR0194

[www.procaredmedfurn.com](http://www.procaredmedfurn.com)

### KEY FEATURE

**STRONG CONSTRUCTION:**

THE FOOT SUCTION MACHINE COMES WITH A POWDER-COATED MILD STEEL BODY, WHICH IS STRONG, DURABLE, AND PROVIDES LONG-LASTING PERFORMANCE.

**SUPERIOR SUCTION JARS:**

THE FOOT SUCTION APPARATUS HAS A PAIR OF 600 ML PVC JARS, FOR A PROLONGED OPERATION AND SUFFICIENT CAPACITY.

**STRONG VACUUM PRESSURE:**

THE DEVICE HAS AN EFFECTIVE AND STRONG SUCTION CAPACITY WITH A VACUUM PRESSURE OF 50 MM HG, TO ENSURE EFFECTIVE FLUID EVACUATION.

**WEIGHT OF THE MACHINE:**

THIS FOOT OPERATED SUCTION MACHINE WEIGHS JUST 4.5 KG, WHICH MAKES THIS MACHINE LIGHTWEIGHT AND PORTABLE AND EASY TO CARRY ANYWHERE IN EMERGENCIES

**FOOT OPERATED:**

THE BEST SOLUTION FOR REMOTE AREAS WHERE LIGHT IS UNAVAILABLE, THIS MACHINE COMES WITH A FOOT-OPERATED MECHANISM

**COMPACT DESIGN:**

THE MACHINE HAS A COMPACT DESIGN MEASURING ONLY 25 X 22 X 15 CM, THE MACHINE IS SMALL, PORTABLE, AND NEEDS LOW SPACE.




Jar



Vacuum Pump




CE



# Procure

Address: Gate no 197, Dehu Alandi road, Jadhavwadi, Chikali Pune, Maharashtra-411062  
Mob no: +91 8377832832  
Web Site : [www.procaredmedfurn.com](http://www.procaredmedfurn.com)  
GST NO: 27ABCFP8176F1ZY



Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procure