

Accompanying the Recovery Journey

Create a Safer Environment

4 FOLD CANVAS STRETCHER

SKU: PRO092

www.procaremedfurn.com

- Measures when unfolded: L 221 x W 53 x H 15 cm
- Folded dimensions: L 61 x W 30 x H 12 cm

KEY FEATURE

- THIS FOLDABLE CANVAS STRETCHER COMES WITH A 4 FOLD MECHANISM WHICH MAKES IT COMPACT AND PORTABLE TO CARRY IN ANY SITUATION TO RESCUE PATIENTS FROM EMERGENCY AREAS
- **STURDY CONSTRUCTION:**
- THE LIGHTWEIGHT AND STRONG FRAME IS MADE FROM PREMIUM MILD STEEL COATED WITH ZINC. MAKE THIS CANVAS STRETCHER MORE DURABLE, AND CORROSION RESISTANT.
- **VINYL-COATED SYNTHETIC FABRIC:**
- THE SURFACE OF THE FOLDABLE CANVAS STRETCHER IS MADE FROM HYGIENIC VINYL-COATED SYNTHETIC FIBER TO PROVIDE THE BEST COMFORT AND SAFETY TO THE PATIENT
- **HOLDING STRAPS:**
- THE NYLON STRAPS HELP TO HOLD THE PATIENT ON THE STRETCHER PREVENT ACCIDENTAL FALLS AND ENSURE PATIENT SAFETY IN TRANSPORTATION.




Stand



Grab
Handle




CE



ProCare

Address: Gate no 197, Dehu Alandi road,
Jadhavwadi, Chikali Pune,
Maharashtra-411062
Mob no: +91 8377832832
Web Site : www.procaremedfurn.com
GST NO: 27ABCFP8176F1ZY



Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

ProCare