

Accompanying the Recovery Journey

Create a Safer Environment

KICK BOWL

SKU: PRO129


www.procaremedfurn.com

- Diameter: **350 mm**
- Height: **280 mm**
- Material – **Stainless steel**
- Wheels – **Comes with the Three 2 inch Wheels**

KEY FEATURE

- KICK BOWL FEATURES
- THE KICK BOWLS ARE DESIGNED TO MAINTAIN HYGIENE DURING MEDICAL OR SURGICAL PROCEDURES. IT COMES WITH 350 MM DIAMETER AND 280 MM HEIGHT WHICH OFFER ENOUGH SPACE AND CAPACITY TO STORE SURGICAL WASTE, HEIGHT MAKE IS EASY TO ACCESS IN EMERGENCIES. THE STAINLESS STEEL CONSTRUCTION OF THE KICK BUCKET MAKES IT CORROSION-RESISTANT, DURABLE TO USE, AND EASY TO CLEAN.
- THE KICK BOWL IS MOUNTED ON THREE 2-INCH WHEELS WHICH HELP TO MOVE THE BOWL EASILY ON DIFFERENT SURFACES, THIS FUNCTION OF THE SURGICAL KICK BUCKET MAKES IT MORE USABLE AND FUNCTIONAL TO USE IN A HOSPITAL TO COLLECT SURGICAL WASTE.
-

stainless steel



Wheel



CE

Procare

Address: Gate no 197, Dehu Alandi road,
Jadhavwadi, Chikali Pune,
Maharashtra-411062
Mob no: +91 8377832832
Web Site : www.procaremedfurn.com
GST NO: 27ABCFP8176F1ZY

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare