



Accompanying the Recovery Journey

MATERNITY PLAIN BED WITH POWDER COATED CRADLE AND BOWS

SKU: PRO019

www.procaremedfurn.com

- Length: 1830 mm
- Width: 830 mm
- Height: 560 mm

KEY FEATURE

Procare's Plain maternity bed comes with powder-coated bows, multi-layer coating enhances the product's durability, and look and offers safety to the patient.

This plain maternity bed offers premium construction, the bed is Pretreated and epoxy powder-coated which offers excellent durability

This plain maternity bed comes with an inbuilt PC baby cradle attachment to provide the best comfort and safety to the newborn baby, the cradle is placed in front of the mother which helps mother to observe her baby while resting

The Two I.V. Rod positions offer multiple options to attach I.V. Rod as per requirements. it also improves access to the I.V. Rod in emergencies.

The provision of the Urine Bag holder improves the health staff's convenience and patient comfort. In addition improves hygienic environment of the ward

Create a Safer Environment

Powder Coated



Bows



Cradle



Procare

Address: Gate no 197, Dehu Alandi road,
Jadhavwadi, Chikali Pune,
Maharashtra-411062
Mob no: +91 8377832832
Web Site : www.procaremedfurn.com
GST NO: 27ABCFP8176F1ZY

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare