

# Accompanying the Recovery Journey

Create a Safer Environment

## WALL MOUNTED MONITOR STAND

SKU: PRO105

[www.procaremedfurn.com](http://www.procaremedfurn.com)

- Length: 500 mm
- Width: 400 mm
- Height: 750 mm

## KEY FEATURE

- THIS WALL-MOUNTED MONITOR IS LIGHTWEIGHT AND IS MADE FROM ALUMINUM ALLOY TO OFFER THE BEST STRENGTH AND DURABILITY.

THIS PARA MONITOR WALL-MOUNTED STAND OFFERS 360° HORIZONTAL RATIIONS TO PROVIDE AN ACCURATE MONITOR VIEW CONVENIENCE DURING TREATMENTS AND OBSERVATIONS.

THE MONITOR STAND OFFERS A 15° VERTICAL TILT FOR EASY UP AND DOWN TO PROVIDE OPTIMAL ANGLE AND OFFERS COMFORT TO VIEWERS DURING OBSERVATION BY REDUCING PRESSURE ON THE NECK AND EYES.

THE PLAITS OF THE MONITOR COME WITH THE SELF-LOCKING RAILS SYSTEM. WHICH HELPS IN EASY INSTALLATION AND ATTACHES THE MONITOR PERFECTLY. THIS SHELF-LOCKING MECHANISM HELPS TO KEEP THE MONITOR STEADY AND SAFE IN ITS PLACE

A 30 CM WALL CHANNEL SUPPORT IS PROVIDED WITH THE STAND TO KEEP THE MONITOR STABLE, IT ENSURES THE STABILITY OF THE MONITOR WHEN THE MONITOR AND ACCESSORIES ARE ATTACHED TO THE STAND. THE WEIGHT CAPACITY IS 20 KG

THE IS IMPORTANT TO MOUNT THE MONITOR AND ACCESSORIES AT THE ACCURATE HEIGHT AND POSITION TO OBSERVE THE MONITOR FROM TIME TO TIME WITHOUT ANY DIFFICULTIES. THE MONITOR COMES WITH A 300 MM HEIGHT


A PREMIUM METAL STORAGE BASKET IS PROVIDED WITH A WALL-MOUNTED MONITOR STAND TO KEEP EQUIPMENT LIKE BP CUFFS, EMERGENCY EQUIPMENT, ETC, AND HELPS TO KEEP THINGS ORGANIZED AND EASILY ACCESSIBLE.

Aluminum Alloy

Wall  
Mount



Basket



CE

# Procare

Address: Gate no 197, Dehu Alandi road,  
Jadhavwadi, Chikali Pune,  
Maharashtra-411062

Mob no: +91 8377832832

Web Site : [www.procaremedfurn.com](http://www.procaremedfurn.com)

GST NO: 27ABCFP8176F1ZY

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare