

Accompanying the Recovery Journey

Create a Safer Environment

PLAIN HOSPITAL BED

SKU: PR0015

www.procaremedfurn.com

Overall Size –

- Length: **1830 mm,**
- Width: **830 mm**
- Height: **560 mm**

THE PLAIN BED IS MADE FROM EPOXY POWDER COATED, PRETREATED CRCA STEEL, THIS PLAIN HOSPITAL BED IS CORROSION RESISTANT AND PROVIDES A LONG LIFE.

THIS PLAIN BED COMES WITH TWO I.V. ROD PROVISIONS WHICH PROVIDE MORE FLEXIBILITY, AND MAKE THIS BED MORE CONVENIENT TO USE IN DAILY ROUTINE.

THE BUILT-IN URINE BAG HOLDER ADDS MORE FUNCTIONALITY, IT HELPS TO MAINTAIN HYGIENE IN WARDS IN DAILY PROCEDURES AND PATIENT CARE.

THESE PLAIN HOSPITAL BED SIZES ARE CUSTOMIZABLE AS PER THE SPECIFIC NEEDS OF THE USERS


THE LENGTH OF THE TABLE IS 1880MM, WIDTH 830 MM AND HEIGHT 570MM WHICH PROVIDE BETTER SPACE AND COMFORT TO THE PATENT



Bows




CE



Procare

Address: Gate no 197, Dehu Alandi road,
Jadhavwadi, Chikali Pune,
Maharashtra-411062
Mob no: +91 8377832832
Web Site : www.procaremedfurn.com
GST NO: 27ABCFP8176F1ZY



Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare