



Procure

Accompanying the Recovery Journey

Create a Safer Environment

PORTABLE SUCTION
MACHINE
SKU: PRO201
www.procuremedfurn.com

- Power Input:** 90 VA
- Energy Source:** DC 12 V, 44 VA
- Maximum Negative Pressure:** ≥ 0.075 MPa or less
- Pumping Rate:** 7E-A ≥ 17 L/min
- Noise level of operation:** < 65 dB (A)
- Jar Capacity:** One portion, 1000 mL of reservoir
- Mode of Operation:** Automatic Pump
Piston Compressor in type

KEY FEATURE

STRONG PERFORMANCE:
THE MACHINE COMES WITH A PISTON COMPRESSOR PUMP, WITH A MAXIMUM NEGATIVE PRESSURE VALUE OF ≥ 0.075 MPa TO PROVIDE AN EFFECTIVE SUCTION WITH A STRONG PUMPING RATE OF ≥ 17 L/MIN

SILENT OPERATION:
TO MAINTAIN THE SILENCE OF THE HOSPITAL AND TO CONTROL NOISE LEVEL FOR BETTER PATIENT COMFORT AND CONCENTRATION, WITH AN OPERATING NOISE LEVEL OF ≤ 65 DB (A)

HIGH CAPACITY:
THE MACHINE HAS A SINGLE POLYCARBONATE JAR WITH A CAPACITY OF 1000 ML THAT OFFERS MORE STORAGE AND LONG USE.

AUTOMATIC OPERATION:
THIS PORTABLE SUCTION APPARATUS COMES WITH AN AUTOMATIC MODE, ITS AUTOMATIC OPERATION MODE ENABLES SMOOTH AND HASSLE-FREE OPERATION

ACCURATE CONTROL:
THE VACUUM PRESSURE METER IS PROVIDED, FOR AN ACCURATE AND CONTROLLED PROCEDURE.

FREE OIL PUMP:
THE MACHINE HAS AN OIL-FREE AND MAINTENANCE-FREE SUCTION PUMP PROVISION FOR BETTER AND PROLONGED WORK EFFICIENCY.

FLEXIBLE POWER OPTIONS:
THE MACHINE IS USABLE IN DIFFERENT AREAS, THIS PORTABLE MACHINE IS POWERED BY A DC 12 V SUPPLY WITH AN INPUT POWER OF 90 VA AND POWER SUPPLY OF 44 VA.

Jar

Adjustable knob



CE

Procure

Address: Gate no 197, Dehu Alandi road,
Jadhavwadi, Chikali Pune,
Maharashtra-411062
Mob no: +91 8377832832
Web Site : www.procuremedfurn.com
GST NO: 27ABCFP8176F1ZY

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procure