



# Accompanying the Recovery Journey

## MATERNITY SEMI FOWLER BED WITH SS CRADLE AND BOWS

SKU: PRO020

[www.procaremedfurn.com](http://www.procaremedfurn.com)

- Length: 1880 mm
- Width: 890 mm
- Height: 560 mm

### KEY FEATURE

You can easily adjust the backrest of the bed for better positioning. The backrest of the semi-fowler maternity bed can rise to 70° to provide the best support to the mother.

Procare's semi-fowler bed comes with a stainless steel baby cradle and bows to provide comfort to mom and baby in any maternity ward. Its stainless steel material makes the cradle more durable, and corrosion-resistant

The bed is made of premium material and has an epoxy powder-coated surface which offers excellent durability and a shiny look. The bed comes with a Stainless steel head and foot board which make it more strong and flexible.

This semi-fowler bed is equipped with a stainless steel baby cradle, the cradle is attached to the foot board and provides the best safety and security to the baby.

The bed is designed to provide two I.V poles at different locations which is easy to access.



Crank



Bows



Cradel



Create a Safer Environment

# Procare

Address: Gate no 197, Dehu Alandi road,  
Jadhavwadi, Chikali Pune,  
Maharashtra-411062  
Mob no: +91 8377832832  
Web Site : [www.procaremedfurn.com](http://www.procaremedfurn.com)  
GST NO: 27ABCFP8176F1ZY

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare