



Create a Safer Environment

SEMI ELECTRIC OPERATING TABLE

SKU: PROOT003
www.procaremedfurn.com

- **Top Dimension-** Length 1980 mm X Width 530 mm
- **Height Adjustment-** 730 mm to 1030 mm
- **Reverse Trendelenburg-** 20° / 25°
- **Lateral Tilt-** 20° / 20°
- **Kidney Elevator-** 130 mm
- **Back Rest (Up and Down)-** 80° / 35°
- **Leg Rest- (Up / Down)-** 10° / 90°
- **Head Rest (Up and Down)-** 45° / 90°
- **Patient Weight Capacity-** 200 Kgs

KEY FEATURE

During procedures, improved imaging and diagnostics are made possible by specially engineered eccentric base, and trouble-free motions of the C-arm.

Our table is outfitted with strong and long-lasting pillars and Premium stainless steel fitting that provides stability and long life to the product

Our semi-electric bed allows interchanging its head and leg portions for better comfort and flexibility of the OT table

The kidney bridge plays an essential role in different surgical procedures, this OT Table comes with integrated kidney bridge option.

This semi-electric Operating table comes with a removable pelvis, leg, and head section which helps healthcare professionals adjust OT table as per the requirement and more comfort while performing the surgeries.

To maintain the Hygiene of the bed and OT room, the table is built using a stainless steel and washable cover.

Adjustable
Head Rest

Shoulder
Support

Waist
Support

Height
Adjustable

Premium Quality
Cushion

Leg Support



CE

ProCare

Address: Gate no 197, Dehu Alandi road,
Jadhavwadi, Chikali Pune,
Maharashtra-411062
Mob no: +91 8377832832
Web Site : www.procaremedfurn.com
GST NO: 27ABCFP8176F1ZY

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

ProCare