



Accompanying the Recovery Journey

Create a Safer Environment

ROUND FRAME SMALL OXYGEN TROLLEY

SKU: PRO139

www.procaredmedfurn.com

- Length: 200 mm
- Width: 710 mm
- Height: 1170 mm

KEY FEATURE

PREMIUM QUALITY CONSTRUCTION:

THIS SMALL ROUND-FRAME OXYGEN CYLINDER TROLLEY IS PRETREATED AND EPOXY POWDER COATED AND MADE FROM PREMIUM CRCA STAINLESS STEEL, AND OFFERS CORROSION-RESISTANT, SHINY, AND STRONG BUILD QUALITY FOR REGULAR USE IN HOSPITALS AND CLINICS

SMOOTH WHEELS:

THIS ROUND FRAME OXYGEN CYLINDER TROLLEY IS EQUIPPED WITH 100-DIAMETER WHEELS, ITS FOUR WHEELBASE MAKES IT MORE SMOOTH AND STABLE ON ANY FLOOR

UNIQUE DESIGN:

THIS OXYGEN CART HAS A ROUND FRAME DESIGN WHICH HELPS TO HOLD THE TROLLEY FROM ALL SIDES, THE SCREWS ARE PROVIDED TO THE OXYGEN CART SO YOU CAN TIGHT THE SCREW TO HOLD THE OXYGEN CYLINDER TIGHTLY.



Oxygen
Cylinder
Support



Wheel



Powder Coated



CE

Procare

Address: Gate no 197, Dehu Alandi road,
Jadhavwadi, Chikali Pune,
Maharashtra-411062
Mob no: +91 8377832832
Web Site : www.procaredmedfurn.com
GST NO: 27ABCFP8176F1ZY

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare