

Accompanying the Recovery Journey

Create a Safer Environment

SMALL OXYGEN TROLLEY

SKU: PRO138

www.procaredmedfurn.com

- Length: 200 mm
- Width: 200 mm
- Height: 1170 mm

KEY FEATURE

WHEELS:

THIS SMALL OXYGEN CYLINDER TROLLEY IS EQUIPPED WITH 100 MM DIAMETER WHEELS WHICH OFFER THE BEST AND EASIEST TRANSPORTATION ON ANY SURFACE.

CYLINDER HOLDERS:

THE CYLINDER HOLDER IS PROVIDED TO HOLD THE CYLINDER ON THE TROLLEY DURING TRANSPORTATION. THIS AVOIDS ACCIDENTAL FALL OF THE CYLINDER FROM TROLLEY

DURABLE CONSTRUCTION:

THE TROLLEY IS MADE FROM PREMIUM CRCA STEEL AND LAYERED WITH AN EPOXY POWDER COATING TO MAKE IT MORE DURABLE, CORROSION-RESISTANT, AND EASY TO CLEAN.



Grab
handle



Safety
Support



Wheel




Powder
Coated



CE

Procare

Address: Gate no 197, Dehu Alandi road,
Jadhavwadi, Chikali Pune,
Maharashtra-411062
Mob no: +91 8377832832
Web Site : www.procaredmedfurn.com
GST NO: 27ABCFP8176F1ZY



Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare