

Accompanying the Recovery Journey

Create a Safer Environment

STRETCHER ON TROLLEY DELUXE SKU: PR0085 www.procaremedfurn.com

- Length: **2100 mm**
- Width: **710 mm**
- Height: **810 mm**
- Body material: The stretcher on a trolley is fully made from Stainless steel to provide hygienic and long-lasting properties
- Side Railings: Complete Stainless Steel sideways swing rails to provide extra safety and simple patient access.
- Wheels: equipped with high-quality Caster wheels come with a **125 mm** diameter
- Shelf: At the bottom of the stretcher trolley one-half shelf is provided for extra storage.
- The oxygen cylinder cage
- Recliner System: This stretcher on trolley comes with a backrest to provide more support and comfort to the patient
- I.V. Pole

KEY FEATURE

- **VERSATILE USAGE:**
- **THIS STRETCHER IS PERFECT TO USE IN DIFFERENT EMERGENCY DEPARTMENTS LIKE AMBULANCES, HOSPITALS, ETC.**
- **NOTE: THE ACTUAL PRODUCT MAY VARY THE IMAGE IS FOR REFERENCE ONLY.**

HYGIENE AND STRONG CONSTRUCTION:
THIS STRETCHER IS MADE FROM PREMIUM STAINLESS STEEL WHICH OFFERS THE BEST DURABILITY, RESIST CORROSION, AND INCREASES THE LIFE OF THE PRODUCT

IMPROVED MOBILITY:
THIS DELUXE STRETCHER OFFERS HIGH-QUALITY CASTER WHEELS TO MOVE THE TROLLEY EASILY AND KEEP IT STABLE.


SAFETY AND CONVENIENCE:
STAINLESS STEEL SIDE RAILING HELPS TO PROVIDE ADDITIONAL SECURITY TO THE PATIENT DURING TRANSPORT



CE

Procare

Address: Gate no 197, Dehu Alandi road,
Jadhavwadi, Chikali Pune,
Maharashtra-411062
Mob no: +91 8377832832
Web Site : www.procaremedfurn.com
GST NO: 27ABCFP8176F1ZY



Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare