

# Accompanying the Recovery Journey

Create a Safer Environment

## SS STRETCHER ON TROLLEY GENERAL

SKU: PRO088

[www.procaremedfurn.com](http://www.procaremedfurn.com)

- Width: 610 mm
- Length: 2100 mm
- Height: 810 mm

### KEY FEATURE

- ITS HIGH-QUALITY STAINLESS STEEL CONSTRUCTION MAKES IT CORROSION-RESISTANT AND MORE DURABLE TO USE.
- DESIGN FEATURES:
- THE TOP OF THE STRETCHER IS MADE FROM STAINLESS STEEL AND IT IS DETACHABLE AND EASY TO CLEAN TO MAINTAIN HYGIENE.
- MOBILITY & STABILITY:
- FOR BETTER MOVEMENT OF THE TROLLEY AND TO PROVIDE BETTER STABILITY ON THE FLOOR 100 MM DIAGONAL WHEELS ARE PROVIDED WITH THE STAINLESS STEEL STRETCHER



Removable  
Stretcher




Grab  
Handle



Wheel




CE



# Procare

Address: Gate no 197, Dehu Alandi road,  
Jadhavwadi, Chikali Pune,  
Maharashtra-411062  
Mob no: +91 8377832832  
Web Site : [www.procaremedfurn.com](http://www.procaremedfurn.com)  
GST NO: 27ABCFP8176F1ZY



Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare