



# Accompanying the Recovery Journey

## YUWELL 7A-23D ELECTRIC SUCTION MACHINE

SKU: PRO196

[www.procaremedfurn.com](http://www.procaremedfurn.com)

- Machine Needs AC power input 180 VA at 220 V at 50 Hz for efficient working.
- The maximum pumping capacity of this suction apparatus is 20 LPM 90 kPa is the maximum negative pressure.
- The noise level of this suction machine is < 65 dB
- Two Jars are provided with 2.5 liter capacity.

## KEY FEATURE

**OIL-FREE, MAINTENANCE-FREE SUCTION PUMP:**  
THIS SUCTION 7A-23D SUCTION APPARATUS PROVIDES OIL FREE PUMP FOR A LONG PERFORMANCE.

**20 LPM CAPABILITY:**  
THESE SUCTION MACHINES HAVE A POWERFUL 20 LITERS PER MINUTE PUMPING CAPACITY TO PROVIDE A SMOOTH AND EFFICIENT SUCTION.

**ADJUSTABLE NEGATIVE PRESSURE:**  
YOU CAN EASILY MANAGE THE MAXIMUM NEGATIVE PRESSURE WITH THE REGULATOR PROVIDED ON THE FRONT PANEL, I WILL GIVE FULL CONTROL TO MANGE AND REGULATE PRESSURE AS PER THE SURGICAL NEED

**STURDY AND STYLISH DESIGN:**  
OUR 7A-23D SUCTION MACHINES ARE MADE OF PREMIUM ABS PLASTIC AND POWDER COATED ALUMINUM, IT HAS FOLDABLE HANDLES THAT MAKE THEM MORE COMPACT AND SIMPLE TO TRANSPORT

**GLASS JARS:**  
THE COMPLETE READY-TO-USE SETUP IS PROVIDED WITH THE MACHINE ITS ACCESSORY KIT CONTAINS TWO 2.5 LITER CAPACITY GLASS JARS, A TUBING KIT, AND AN EXTRA BACTERIAL FILTER.

Remote



Jar



Grab Handle



Wheel



Create a Safer Environment

CE

# Procare

Address: Gate no 197, Dehu Alandi road,  
Jadhavwadi, Chikali Pune,  
Maharashtra-411062  
Mob no: +91 8377832832  
Web Site : [www.procaremedfurn.com](http://www.procaremedfurn.com)  
GST NO: 27ABCFP8176F1ZY

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare