

Accompanying the Recovery Journey

Create a Safer Environment

YUWELL 7A-23A SUCTION MACHINE

SKU: PRO196

www.procaremedfurn.com

- The required power supply for the Suction machine is AC 220 V at 50Hz, 180 VA
- The Pumping capacity of the suction apparatus is 20 LPM
- You can increase the negative pressure of this machine up to 90 kPa
- Working Noise: Less than 60 dB
- Two Polycarbonate Jar with a Capacity of: 2.5 L x 2

KEY FEATURE

OIL-FREE SUCTION PUMP:

OUR SUCTION MACHINE COMES WITH ADVANCED OIL-FREE PUMPS THAT PROVIDE A BETTER PERFORMANCE AND DELIVER RESULTS.

20 LPM CAPABILITY:

YUWELL 7A-23A SUCTION HAS A 20 LITER PER MINUTE PUMPING CAPACITY OR A QUICK AND EFFECTIVE SUCTION.


ADJUSTABLE NEGATIVE PRESSURE:

ON THE FRONT PANEL OF THE SUCTION MACHINE REGULATORS ARE PROVIDED TO CONTROL THE NEGATIVE PRESSURE OF THE MACHINE

ROBUST AND STYLISH DESIGN:

OUR SUCTION APPARATUS COMES WITH AN ABS BODY AND FOLDING HANDLES, THIS PREMIUM ABS MATERIAL MAKES THIS MACHINE MORE DURABLE AND ROBUST.

Jar



Wheel



Procare

Address: Gate no 197, Dehu Alandi road,
Jadhavwadi, Chikali Pune,
Maharashtra-411062
Mob no: +91 8377832832
Web Site : www.procaremedfurn.com
GST NO: 27ABCFP8176F1ZY

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

Procare